

Summary Report of Clinical Study CELLRETURN NECKLE RAY

<p>Results of Study</p>	<p>This clinical study was processed with females who are aged between 29 and 62 and suitable to subject selection criteria for the test. Subjects used test product "CELLRETURN NECKLE RAY and Therapy Advance Lifting Neck Cream" for 4 weeks. Effects of thick neck wrinkle, fine neck wrinkle, depth of neck wrinkle, sagging neck lifting, skin elasticity, resilience, skin density and skin texture improvement was assessed.</p> <p>1) The total 22 subjects who finished this clinical study were females. The average age was 53.409 years old. Selected subjects did not show any adverse events or any medical or drug history that might affect the study.</p> <p>2) As a results of checking thick neck wrinkle, the value significantly decreased($p<0.05$) after 2 and 4 weeks of product use, compared to before use.</p> <p>3) As a results of checking thick neck wrinkle, the value significantly decreased($p<0.05$) after 2 and 4 weeks of product use, compared to before use.</p> <p>4) As a results of checking thick neck wrinkle, the value significantly decreased($p<0.05$) after 2 and 4 weeks of product use, compared to before use.</p> <p>5) As a results of checking thick neck wrinkle, the value significantly decreased($p<0.05$) after 2 and 4 weeks of product use, compared to before use.</p> <p>6) As a results of checking thick neck wrinkle, the value significantly decreased($p<0.05$) after 2 and 4 weeks of product use, compared to before use.</p> <p>7) As a results of checking thick neck wrinkle, the value significantly decreased($p<0.05$) after 2 and 4 weeks of product use, compared to before use.</p> <p>8) As a results of checking thick neck wrinkle, the value significantly decreased($p<0.05$) after 2 and 4 weeks of product use, compared to before use.</p> <p>9) As a results of checking thick neck wrinkle, the value significantly decreased($p<0.05$) after 2 and 4 weeks of product use, compared to before use.</p> <p>10) As a result of efficacy survey, the 100.000% of subjects have positively assessed as effects of effects of thick neck wrinkle, fine neck wrinkle, depth of neck wrinkle, sagging neck lifting, skin elasticity, resilience, skin density and skintexture improvement.</p> <p>11) There is no report of special adverse event to the skin while the test participant used the test product, and there is no abnormal observation upon the physical examination by a dermatologist.</p> <p>Therefore, it is deemed that the use of "CELLRETURN NECKLE RAY" and "Therapy Advance Lifting Neck Cream" for 4 weeks will help improve thick neck wrinkle, fine neck wrinkle, depth of neck wrinkle, sagging neck lifting, skin elasticity, resilience, skin density and skin texture</p>
-------------------------	--